CARE FLOW Foley-Katheter

D Foley-Katheter

Produktbeschreibung

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GB Foley-Catheter

Device description

Flexible shaft made of multiple-cleansed latex with siliconised or silicone elastomer coated surface

Indications

- acute an chronic urinary disorders
- intra- and postoperative urinary drainage
- urinary drainage for measurement of urine excretion e.g. during intensive care and infusion therapy

Directions for use

- The device is sterile as long as the packaging using aseptic technique.
- Remove the catheter from its sterile packaging using aseptic technique..
- Using sterile technique as indicated in published literature, insert the catheter fully into the urethra ensuring that the balloon is beyond the bladder neck.
- Inflate the balloon with sterile water. Use sterile water only.
- The capacity of the balloon is indicated on the funnel of the catheter or on the label of the packaging. Do not overinflate.
- To deflate the balloon prior to removal aspirate the syringe gently to remove the inflation fluid. Do not use excessive aspiration on the syringe during deflation as this may cause a vacuum collapse of the inflation lumen which may impair normal drainage.
- Maximum recommended indwelling times will vary due to patient specific factors.

Potential complications

If difficulty is encountered aspirating the balloon with the syringe, a rare and infrequently reported event, the limb of the catheter with the valve should be cut with sharp scissors at the bifurcation, or the balloon ruptured according to established procedures reported in medical literature. Should it be necessary to rupture the balloon, care must be taken to remove all fragments from patient. Irritation of the urethral mucosa, blockage of the catheter due to encrustation and catheter induced infections are documented complications with some catheter materials and patients.
The patient should be routinely monitored in accordance with accepted procedures.
The siliconised catheter should be changed after an appropriate time usually not more than 7 days or specified by a medical practitioner.
Usage insitu for silicone elastomer coated catheters should not be more than 4 weeks or specified by a medical practitioner.

Warnings and Precautions

For urological use only.
Device contains latex which may elicit allergic responses in latex sensitized individuals.
Latex Urethral Catheter is packaged for single use.
Do not reuse, reprocess. Reuse, reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient.
Reuse, reprocessing may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.
Do not use lubricants or preparations with a petroleum base on devices which contain latex.
Always inflate the balloon with sterile water.
Never clamp the catheter.
Never use a syringe to puncture the catheter shaft for urine sampling. There is a danger of needle-stick injury and the catheter may be compromised.

Excessive aspiration can collapse inflation lumen, preventing balloon deflation.
Visually inspect the package for any breach of the packaging integrity.
For device complaints, please contact your supplier. Please also report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority at the user's place of residence.

Paediatric Catheters

Paediatric catheters include a stylet / stillette to facilitate insertion. Make sure stillette end is properly positioned in catheter tip and not protruding from the catheter eye during insertion.
Remove the stillette once the catheter is positioned correctly before inflating the balloon.

Storage

These devices should be stored in their original box, preferable away from direct and indirect sources of light and heat.
After use, this device may be a potential biohazard. Handle and dispose off in accordance with accepted medical practice and applicable laws and regulations.

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